

OptoTech's 4K Technology builds on the OptoDesign Pro (ODP) platform, using an advanced, **high-resolution design optimization process**. By dramatically increasing the number of iterations in the optimization algorithm, 4K Technology delivers **exceptionally smooth power transitions** with ultra-high optical precision. The result for the wearer: brilliant, crisp, and clear vision across a wide field of view.



4K Benefits:

- + *Advanced design optimization method*
- + *More iteration steps for smooth power gradients and very high optical resolution*
- + *Brilliant sharp and clear vision*
- + *Wide optical fields*