

Lens Comparison Chart			DESIGN DETAILS & IDEAL WEARER	COMPARE TO	
ACCULAB LENSES					
PROGRESSIVES	BEST	MONTANEL 	High end, all purpose, automatic corridor · Seasoned PAL wearer		• Essilor Varilux X Series & W3+ • Shamir (Essilor) Autograph Intelligence 2.0 • Hoyalux iD MyStyle V, iD MyStyle 3 & LifeStyle 4 • Zeiss Individual 3
		IOT CAMBER STEADY PLUS (w/ AR)	High end, all purpose, Poseidon or Neptune Plus AR included, automatic or manual corridor · Seasoned PAL wearer		
		SEIKO SUPERIOR	Premium, all purpose, auto/manual corridor · Experienced PAL wearer		
	BETTER	MONTANEL 	Premium, all purpose, automatic corridor · Experienced PAL wearer, also good for patients transitioning from SV, bifocal, or anti-fatigue lenses		• Essilor Physio & Comfort W2 • Hoyalux iD LifeStyle 3, LifeStyle 4, MySelf, & Array 2 • Shamir (Essilor) Autograph II & Autograph IIID • Zeiss SmartLife Individual & SmartLife 2/3
		IOT ENDLESS STEADY	"Better", all purpose, auto or manual corridor · Regular PAL wearer		
		SEIKO SURMOUNT · SUPERCEDE II			
	GOOD	MONTANEL 	Entry-level/Economical PAL, automatic corridor · For patients transitioning from SV, bifocal, or anti-fatigue lenses or regular PAL wearers		• Essilor Varilux Comfort Max • Hoyalux Balansis • Shamir (Essilor) HD Experience, Element & Genesis HD • Zeiss SmartLife
		IOT ESSENTIAL STEADY	Entry-level/Economical PAL, auto or manual corridor · For patients transitioning from SV, bifocal, or anti-fatigue lenses		
		SEIKO SUPERNAL			
	NEAR	MONTANEL 	Indoor lens for office/reading, optimized for Near vision, automatic corridor · Patients looking for a PAL w/ focus on the near zone		• Essilor Varilux X Series & Comfort Short • Shamir (Essilor) Autograph II+
		SEIKO 	Zero drop, all purpose, auto/manual corridor · Good for plus patients		
FAR	MONTANEL 	Outdoor lens for sports/driving, optimized for Distance vision, auto/manual corridor · Patients looking for a PAL w/ focus in the far zone		• Hoya iD LifeStyle 3 Outdoor & MyStyle 3 Adventure	
LIFESTYLE	OFFICE	MONTANEL 	Computer/Office lens, optimized for Near & Intermediate, automatic corridor · Office workers		• Essilor Computer • Shamir (Essilor) Workspace, WorkStyle 3, InTouch, Computer, SmartOffice, & Autograph II Office • Hoya Tact BKS • Zeiss OfficeLens Series
		IOT ENDLESS OFFICE	Computer/Office lens. Comes in 3 levels based on the distance portion: 2ft - 4ft, 2ft - 6.5ft, & 2ft - 13ft. · Office workers		
	ANTI-FATIGUE	MONTANEL 	Single Vision lens with add boost in reading area, add blue block to match Eyezen · Available in: 0.40, 0.60, 0.85, 1.10, & 1.20 diopters		• Essilor Eyezen+ Series • Shamir (Essilor) Relax • Hoya Sync III • Zeiss SmartLife Digital
		IOT ENDLESS ANTI-FATIGUE	Single Vision lens with add boost in reading area, add blue block to match Eyezen · Available in: 0.0, 0.50, 0.75, 1.00, & 1.25 diopters		