

ACCULAB OPTICAL

ANTI-FATIGUE LENS DISPENSER GUIDE:

AccuLab offers **10 different anti-fatigue lens designs** for your patients. Use the chart below to determine which design is right for them. The difference between the designs is the amount of accommodative relief the patient will receive in the lens:

3 Different Lens Design Options:
(Compare to Eyezen)*

- IOT ENDLESS SINGLE VISION
- IOT ENDLESS ANTI-FATIGUE
- MONTANEL DART FATIGUE SHIELD

*Pair with Trident AR or AccuBlu lens material for blue light protection.

0.00 diopters	0.40 diopters	0.50 diopters	0.60 diopters	0.75 diopters	0.85 diopters	1.00 diopters	1.10 diopters	1.20 diopters	1.25 diopters
IOT ENDLESS SINGLE VISION 0.0 diopters	MONTANEL DART FATIGUE SHIELD 1 0.40 diopters	IOT ENDLESS ANTI-FATIGUE 050 0.50 diopters	MONTANEL DART FATIGUE SHIELD 2 0.60 diopters	IOT ENDLESS ANTI-FATIGUE 075 0.75 diopters	MONTANEL DART FATIGUE SHIELD 3 0.85 diopters	IOT ENDLESS ANTI-FATIGUE 100 1.00 diopters	MONTANEL DART FATIGUE SHIELD 4 1.10 diopters	MONTANEL DART FATIGUE SHIELD 5 1.20 diopters	IOT ENDLESS ANTI-FATIGUE 125 1.25 diopters
Compare to (EYEZEN START) & (EYEZEN+ 0)	Compare to (EYEZEN+ 1)		Compare to (EYEZEN+ 2)		Compare to (EYEZEN+ 3)		Compare to (EYEZEN+ 4)		
Patients age 17 & under	Patients age 18 to 34	Patients age 25 to 35	Patients age 35 to 39	Patients age 35 to 40	Patients age 40 to 44	Patients age 45 to 50	Patients age 45 to 60	Patients age 50 to 65	Patients age 50 to 65