



Chemistrie Color lenses filter precise portions of the visible light spectrum to create contrast and color separation, enabling patients with Color Vision Deficiency (CVD) or Color Blindness to discern additional colors or shades.

Available as a

- Chemistrie Clip
- Uncut Lens
- Readymade Frame + Lens

Three strengths available:

- Indoor
- Outdoor Moderate
- Outdoor Severe
- Strengths decided based upon patient's need.
- Indoor: Not as dark as other lenses. Backside AR
- Outdoor Moderate: Sunglass with a little bit of red hue. Backside AR
- ⬆ *Most likely most popular option*
- Outdoor Severe: Sunglass with a darker hue. Backside AR.



- Uncut**
- Thickness: 1.5mm
 - Available in 4 & 6 Base Curves

- Clip**
- Thickness: 1.2mm
 - Available in 2, 4, 6 Base Curves

Readymade Styles available:

Any of three frame styles, available in any of the three strengths

- **CB10 Bright: Black, Tortoise or Blue**



- Plastic Frame

- **CB11 Light: Gold, Silver or Black**



- Stainless Steel Frame

- **CB12 Lumin: Gold, Silver or Black**



- Stainless Steel Frame

- All lenses are made of Poly
- All Frames come with Case, Box, and Cleaning Cloth
- Uncut Lenses also come with Case, Box, and Cleaning Cloth

POP
\$285

Includes:

- Display
- 3 Demos Frames (2 Readymade Styles, One Clip)
- 3 Readymade Styles to sell (One of each frame style/lens strength)

*Standard Assortment Comes with the Display, not a customizable option

