











Transitions® Light Intelligent Lenses™ by Feature

	 Our Best Everyday Lenses	 Our Best Lenses for Extra Protection, Even in the Car	 Our Best Lenses for Glare Reduction Outdoors	 Adaptive Sunglass Lenses for Behind the Wheel
ICONIC LENS COLORS	 Graphite Green, Gray, or Brown	 Graphite Green, Gray, or Brown	 Gray	 From Olive Green to Copper to Brown
NEW COLORS	 Style Colors: Amber, Emerald, Sapphire, and Amethyst	 Style Mirrors: Gold, Silver Shadow, Pink, Red, Green, and Blue		
INDOOR CLARITY	FULLY CLEAR	HINT OF TINT FOR PROTECTION	HINT OF TINT FOR PROTECTION	OLIVE GREEN <i>Not Recommended for Indoor Use</i>
OUTDOOR DARKNESS	DARK	DARKEST <i>even in high temperatures</i>	DARK	DARK
BEHIND-THE-WHEEL ACTIVATION	NO	YES	NO	YES
POLARIZATION What's this?	NO	NO	YES - VARIABLE POLARIZATION	YES - POLARIZED