

Adaptive Office design  
with slight power boost

Office lenses are the ideal solution for today's workplace, giving you the perfect visual boost for everything from reading and computer work to team meetings. The **Javelin** offers **expanded near and intermediate vision zones**, providing exceptional comfort for close-up and mid-range tasks. Its soft design ensures **smooth transitions** between focal distances, making it easy to adapt from one task to another.

This updated design also features a **variable corridor length**, allowing each pair to be customized to fit the wearer's unique needs. Built with OptoTech's **OptoDesign Pro Technology**, along with their **OMD and 4K Technology**, these lenses deliver sharp, crystal-clear vision with outstanding detail and clarity.



Design Details & Specs	
Lifestyle · Near / Intermediate · Office	
Available corridor length	Automatic
Fitting cross to PRP distance	0 mm
Min. fitting height	16 mm

Prescribed		Dynamic Power Office Lens		
Add. Power	-0.75	-1.25	-1.75	-2.25
0.75	Infinity			
1.00	4.00			
1.25	2.00	Infinity		
1.50	1.35	4.00		
1.75	1.00	2.00	Infinity	
2.00	0.80	1.35	4.00	
2.25		1.00	2.00	Infinity
2.50		0.80	1.35	4.00
2.75			1.00	2.00
3.00			0.80	1.35
3.25				1.00
3.50				0.80

Finer add. power gradations possible

Ideal for

Experienced progressive wearers who are looking for superior clear vision while working

Application

Large near and intermediate visual fields for working with computers & digital displays or other activities (e.g. reading, cooking, painting)

Office  
**JAVELIN**



Near ●●●●○ Intermediate ●●●●● Far ●●○○○ Adaptive ●●●●●

Wearer's Benefits:

- + Natural vision and improved ergonomic posture
- + Enhanced near and intermediate visual fields
- + Easier adapting vision between zones
- + Fully customized lens