

High-performance with generous visual fields and great image quality

The **Rondel** is a soft progressive lens design that offers **smooth transitions between near and far vision**. With wide viewing zones and consistently sharp image quality, it's an excellent choice for **everyday wear**. Designed with advanced **4K Technology**, these lenses deliver clear, comfortable vision and make adapting to them quick and easy.



Design Details & Specs	
Better · Tier 4 Category N · All purpose · Experienced wearer	
Available corridor length	Automatic
Fitting cross to PRP distance	4 mm
Min. fitting height	14 mm

**Ideal for**

Experienced wearers who are looking for an all-purpose lens of high-quality

**Application**

Clear and comfortable view for every day

Premium

**RONDEL**

PROGRESSIVE

LIFESTYLE

SPECIAL

Near ●●●○○

Far ●●●○○

Intermediate ●●●●○

Adaptive ●●●●○

**Wearer's Benefits:**

+ Well balanced in all zones

+ Reduced eyestrain

+ Greater comfort

+ Fully customized lens

+ Easy adaptation between zones