

Balanced progressive lens design with optimized visual fields

OptoTech's **OptoDesign Pro (ODP) Technology** powers our **KBar** progressive lens design to deliver clear, comfortable vision with **evenly balanced visual fields**. In this design, all viewing zones are given equal weight, creating a smooth, natural transition into astigmatism areas in both near and distance vision. **ODP Technology** ensures these areas stay balanced while **minimizing distortion**, so wearers enjoy a more relaxed, natural visual experience.



Design Details & Specs

Good - Tier 3 Category F - All purpose - Entry level

Available corridor length	Automatic
Fitting cross to PRP distance	4 mm
Min. fitting height	16 mm

Ideal for

All kind of office work that requires a balanced progressive lens

Application

Optimized and comfortable balanced view for office work

Standard



Near ●●○○○ Far ●●○○○
Intermediate ●●○○○ Adaptive ●●●○○

Wearer's Benefits:

- + *Balanced in all fields of view*
- + *Greater comfort for office activities*
- + *Fully customized lens*
- + *Easy adaptation between zones*

PROGRESSIVE

LIFESTYLE

SPECIAL