

# ATLANTIS™

## SERIES

PREMIUM ANTI-REFLECTIVE COATINGS

One Premium Standard. Four Signature Finishes.



### REDUCES GLARE

For clearer, more comfortable vision



### RESISTS SCRATCHES

Advanced protection for everyday wear



### REPELS WATER & OIL

Easy to clean, stay cleaner longer



### RESISTS SMUDGES

Keeps lenses looking clear and new



### UV PROTECTION

Complete 360° front & back



### BLUE LIGHT PROTECTION

Available with **Trident** coating

## COMPARE OUR COATINGS



### POSEIDON

Green Reflex

#### Best for:

Everyday premium eyewear

- ✓ Water & Oil Repellent
- ✓ UV Protection
- ✓ Scratch Resistant
- ✓ Easy Clean

Comparable to Crizal® Rock



### NEPTUNE

Blue Reflex

#### Best for:

Blue cosmetic reflex

- ✓ Water & Oil Repellent
- ✓ UV Protection
- ✓ Scratch Resistant
- ✓ Easy Clean

Comparable to Crizal® Sapphire



### MIRAGE

Nearly invisible, Subtle Blush Undertones

#### Best for:

Patients wanting the most natural-looking lenses

- ✓ Lowest visible reflex
- ✓ Water & Oil Repellent
- ✓ UV Protection
- ✓ Scratch Resistant
- ✓ Easy Clean

Comparable to Crizal® Natural Look



### TRIDENT

Blue Light Protection, Blue-Violet Reflex

#### Best for:

Digital device users

- ✓ Blue Light Filter
- ✓ Water & Oil Repellent
- ✓ UV Protection
- ✓ Scratch Resistant
- ✓ Easy Clean

Comparable to Crizal® Previncia



## WHY ATLANTIS?



### CLEARER VISION

Advanced AR technology reduces glare for sharper, more comfortable vision.



### BETTER APPEARANCE

Minimizes reflections so others see your eyes, not glare.



### BUILT TO LAST

Superior scratch resistance and easy-clean surface keep lenses looking new longer.



### COMPLETE UV PROTECTION

360° protection on both front and back surfaces.

## ADVANCED TECHNOLOGY



• **EASY-CLEAN TOP COAT**  
Repels water, oil, and dust for easier cleaning

• **HYDROPHOBIC LAYER**  
Repels water and reduces smudges

• **ANTI-REFLECTIVE LAYERS**  
Multiple precision layers reduce reflections for clearer vision

• **ADHESION LAYER**  
Ensures long-lasting durability of the AR stack

• **HARD COAT**  
Provides scratch resistance and durability

• **LENS**  
High-quality optical lens with standard scratch coat

Multiple precision-deposited layers work together to reduce reflections, improve durability, and keep lenses clearer longer.