

## FEATURES



### CLEAR

41 Abbe Value



### THIN

1.60 Index of Refraction



### SAFE

Impact Resistant  
Passes Z87.1-2015 Dropball Test



### EASY TO PROCESS

Cuts like plastic  
Low wear and tear on machines

# Ultra-FLEX<sup>++</sup>



## AVAILABLE IN:

SEMI-FINISHED SV UC  
SEMI-FINISHED SV HC  
SEMI-FINISHED SV BB LENS UV420 HC  
SEMI-FINISHED SV Q-CHANGE HC

MORE AVAILABILITY COMING SOON