

## Balanced progressive lens design with optimized visual fields

OptoTech's **OptoDesign Pro (ODP) Technology** powers our **KBar** progressive lens design to deliver clear, comfortable vision with **evenly balanced visual fields**. In this design, all viewing zones are given equal weight, creating a smooth, natural transition into astigmatism areas in both near and distance vision. **ODP Technology** ensures these areas stay balanced while **minimizing distortion**, so wearers enjoy a more relaxed, natural visual experience.



### Design Details & Specs

**Good · Tier 3 Category F · All purpose · Entry level**

Available corridor length	Automatic
Fitting cross to PRP distance	4 mm
Min. fitting height	16 mm

### Ideal for

All kind of office work that requires a balanced progressive lens

### Application

Optimized and comfortable balanced view for office work

Standard  
**KBAR**



Near ●●○○○ Far ●●○○○  
Intermediate ●●○○○ Adaptive ●●●○○

### Wearer's Benefits:

- + *Balanced in all fields of view*
- + *Greater comfort for office activities*
- + *Fully customized lens*
- + *Easy adaptation between zones*

PROGRESSIVE

LIFESTYLE

SPECIAL

