

FEATURES



CLEAR INDOORS

Over 90% VLT



DARK OUTDOORS

Up to 85% absorption*



FAST CHANGING

Clear to dark in 20 seconds*



UV PROTECTION

Blocks 100% UVA & UVB
Up to 420 nm



BLUE FILTER

Reduces digital eyestrain



HIGH CONTRAST

NeoContrast Technology
Vibrant, vivid, true colors

FUL-PROTECT



AVAILABLE IN:

SEMI FINISHED SV POLY HC - GREY

MORE AVAILABILITY COMING SOON

* Absorption and speed can vary based on lighting conditions, ambient temperature, lens color, and lens material