

Premium

# STILETTO

## High-performance progressive

### for enhanced near vision

The **Stiletto** is a soft, adaptive progressive lens designed to give your patients **a wide, clear view for close-up tasks**. Using OptoTech's **ODP Technology**, the lens layout is optimized so that any distortion is pushed to the outer edges, improving overall visual performance. Paired with **4K Technology**, wearers enjoy high-resolution clarity and an exceptionally **generous near-vision zone**.



Design Details & Specs	
Better · Near - All Purpose Indoor	
Available corridor length	Automatic
Fitting cross to PRP distance	4 mm
Min. fitting height	16 mm

### Wearer's Benefits:

- + Improved near distance focus
- + Reduced eyestrain
- + Greater comfort
- + Fully customized lens
- + Easy adaptation between zones

### Ideal for

Experienced wearers looking for a high-quality lens with improved near vision

### Application

Near distance and indoor activities that requires an open and clear vision (e.g. reading)

