



Digital Single Vision and Near Variable Focus

Progressive Addition Lenses

**ENGRAVING INDEX**

Product	Symbol	VID	MFH	Category	Product	Symbol	VID	Category
Unity Via Office Pro 5	U	Vo5	14	IA	Unity SVx	U	SX	BA
Unity Via Office Pro 10	U	Vo1	14	IA	Unity SVxtra	U	ST	BA
Unity Relieve 50	U	R5	-	BA+TA	Unity SVxtreme	U	SV	BA
Unity Relieve 70	U	R7	-	BA+TA				

- Left lens
- Right lens
- Material Code
- Add Power
- Nose-pad alignment line
- Centration Unity Relieve & Unity Single Vision Lenses
- Unity Via Office Pro

**ENGRAVING INDEX**

Product	Symbol	VID	MFH	Category	Product	Symbol	VID	MFH	Category
Unity Via Elite VR	U	VRe	12	N+TA	Unity Via Plus II	U	VpII	12	0
Unity Via Elite II-Variable	U	VeII	12	N	Unity Via Mobile II	U	VmII	12	0
Unity Via Elite II 12	U	EII2	12	N	Unity Via Wrap II	U	VwII	12	0
Unity Via Elite II 14	U	EII4	14	N	Unity Via II	U	VII	12	F
Unity Via Elite II 16	U	EII6	16	N	Ethos Plus	E	Eps, Ep	14,18	J
Unity Via Elite II 18	U	EII8	18	N	Ethos	E	Et	18	K