

Upgraded premium all-purpose lens with remarkable image stability for excellent vision

The **Arrow** is built for the demanding visual needs of an active lifestyle, delivering wide, crystal-clear vision and exceptional comfort **for everyday activities**. This upgraded design offers major improvements in distance vision and provides an optimized visual experience with **sharper clarity and precision at every range**. With **OMD Technology**, unwanted astigmatism is reduced by up to **40%**, ensuring outstanding visual performance, smooth power transitions between distances, and effortless adaptation.



| Design Details & Specs                          |           |
|---|-----------|
| Best • High End • All purpose • Seasoned wearer |           |
| Available corridor length                       | Automatic |
| Fitting cross to PRP distance                   | 4 mm      |
| Min. fitting height                             | 16 mm     |

**Ideal for**

Wearers looking for a premium quality and enhanced vision at minimal distortion for everyday

**Application**

Enhanced clear and comfortable vision for day to day and for using electronic devices

High End

ARROW

Near

Intermediate

Far

Adaptive

Wearer's Benefits:

+ Superior focus & reduced eyestrain

+ Excellent balanced view

+ Highest comfort

+ Fully customized lens

+ Easy adaptation between zones