

Young style progressives with slight power boost

Not all eyeglass lenses are built for every task. For people who spend long hours reading, working at a desk, or using digital devices, they may need a lens designed to fight eye fatigue.

The **Dart Fatigue Shield** is an **advanced anti-fatigue lens** made for adults who wear single vision lenses and experience tired eyes, especially after extended screen time. It features a **subtle power boost** in the lower portion of the lens that eases eye strain, improves focus, and enhances overall comfort.

With OptoTech's cutting-edge **ODP, 4K, and OMD Technologies**, the Dart provides crisp, high-resolution vision with minimal distortion, **smooth transitions between visual zones**, and effortless adaptation for everyday use.



Design Details & Specs

Lifestyle · Anti-fatigue / Baby Add · Office

Available corridor length	Automatic
Fitting cross to PRP distance	4 mm
Min. fitting height	16 mm

Ideal for

Beginners or non-adapts who are looking for a premium anti-fatigue lens to avoid tired eyes and to improve focus

Application

Comfortable view, especially when using electronic devices like smartphones, tablets, laptops etc.

Lifestyle

DART Fatigue Shield



Near ●●●●○ Far ●●●○○
Intermediate ●●●●○ Adaptive ●●●●●

PROGRESSIVE

LIFESTYLE

SPECIAL

Wearer's Benefits:

- + A slight power boost of a low addition in the lower portion of the lens to reduce eyestrain during close up activities
- + Greater comfort than single vision correction lenses due to the accommodative relief in the near vision



www.acculab.net/montanel

AccuLab Optical