



DESIGN DETAILS & IDEAL WEARER

COMPARE TO

PROGRESSIVES

BEST

MONTANEL **ARROW**

High end, all purpose, automatic corridor · Seasoned PAL wearer

IOT CAMBER STEADY PLUS (w/ AR)

High end, all purpose, **AR included**, automatic or manual corridor · Seasoned PAL wearer

SEIKO SUPERIOR

Premium, all purpose, auto/manual corridor · Experienced PAL wearer

MONTANEL **RONDEL**

Premium, all purpose, automatic corridor · Experienced PAL wearer, also good for patients transitioning from SV, bifocal, or anti-fatigue lenses

IOT ENDLESS STEADY

SEIKO SURMOUNT · SUPERCEDE II

"Better", all purpose, auto or manual corridor · Regular PAL wearer

MONTANEL **KBAR**

Entry-level/Economical PAL, automatic corridor · For patients transitioning from SV, bifocal, or anti-fatigue lenses or regular PAL wearers

IOT ESSENTIAL STEADY

SEIKO SUPERNAL

Entry-level/Economical PAL, auto or manual corridor · For patients transitioning from SV, bifocal, or anti-fatigue lenses

NEAR

MONTANEL **Stiletto**

Indoor lens for office/reading, optimized for Near vision, automatic corridor · Patients looking for a PAL w/ focus on the near zone

SEIKO **SABRE**

Zero drop, all purpose, auto/manual corridor · Good for plus patients

FAR

MONTANEL **LANCÉ**

Outdoor lens for sports/driving, optimized for Distance vision, auto/manual corridor · Patients looking for a PAL w/ focus in the far zone

LIFESTYLE

OFFICE

MONTANEL **JWELIN**

Computer/Office lens, optimized for Near & Intermediate, automatic corridor · Office workers

IOT ENDLESS OFFICE

Computer/Office lens. Comes in 3 levels based on the distance portion: **2ft - 4ft, 2ft - 6.5ft, & 2ft - 13ft.** · Office workers

ANTI-FATIGUE

MONTANEL **dart** Fatigue Shield

Single Vision lens with add boost in reading area, add blue block to match Eyezen · **Available in: 0.40, 0.60, 0.85, 1.10, & 1.20 diopters**

IOT ENDLESS ANTI-FATIGUE

Single Vision lens with add boost in reading area, add blue block to match Eyezen · **Available in: 0.0, 0.50, 0.75, 1.00, & 1.25 diopters**

- Essilor Varilux X Series & W3+
- Shamir (Essilor) Autograph III & Intelligence
- Hoya iD MyStyle 3 & LifeStyle 4

- Essilor Comfort W2
- Hoya Array 2
- Zeiss SmartLife

- Essilor Varilux Comfort Max
- Shamir (Essilor) Element & Genesis HD

- Essilor Varilux X Series & Comfort Short
- Shamir (Essilor) Autograph II+

- Hoya iD LifeStyle 3 Outdoor & MyStyle 3 Adventure

- Essilor Computer
- Shamir (Essilor) Workspace, WorkStyle 3, InTouch, Computer, & Autograph II Office
- Hoya Tact BKS
- Zeiss OfficeLens Series

- Essilor Eyezen+ Series
- Shamir (Essilor) Relax
- Hoya Sync III
- Zeiss SmartLife Digital